

# flexcon<sup>®</sup>

---

C O N T A C T   L E N S E S

**A Guide to Trouble-Free Contact Lens Wear  
(Instructions for Use)**


*Consult an Eye Care Professional for suitability of wear*

[www.oculuslens.com](http://www.oculuslens.com)


Product Type	Contents	Indication of use (Intended Purpose)
Flexcon Daily Wear	One sterile soft contact lens, 38% water, 62% Poly (2-Hydroxyethyl Methacrylate) immersed in 0.9% buffered saline.	Is applied directly to the cornea, and is designed for daily wear for the correction of refractive ametropia (myopia and hyperopia) in aphakic and non-aphakic persons with non-diseased eyes that may exhibit astigmatism up to 1.00 Diopters that does not interfere with visual acuity. Eye Care Professional may prescribe the lens for daily wear on a frequent replacement schedule up to a year.
Flexcon Super-Thin		
Flexcon Blue Tint Super-Thin		
Flexcon Toric Blue Tint	One sterile soft contact lens, 38% water, 62% Poly (2-Hydroxyethyl Methacrylate) immersed in 0.9% buffered saline.	Is applied directly to the cornea, and is designed for daily wear for the correction of refractive ametropia (myopia and hyperopia) in aphakic and non-aphakic persons with non-diseased eyes that may exhibit astigmatism up to 6.00 Diopters. Eye Care Professional may prescribe the lens for daily wear on a frequent replacement schedule up to a year.
Flexcon Extended Wear	One sterile soft contact lens, 55% water, 45% Hydroxyethyl methacrylate-methacrylic acid copolymer immersed in 0.9% buffered saline.	Is applied directly to the cornea, and is designed for daily wear for the correction of refractive ametropia (myopia and hyperopia) in aphakic and non-aphakic persons with non-diseased eyes that may exhibit astigmatism up to 1.00 Diopters that does not interfere with visual acuity. Eye Care Professional may prescribe the lens for daily wear on a frequent replacement schedule up to a year.
Flexcon Prosthetic PRDW	One sterile soft contact lens, 38% water, 62% Poly (2-Hydroxyethyl Methacrylate) immersed in 0.9% buffered saline.	Indicated for concealing unsightly scarring or abnormalities of the iris. The lens may or may not have optical correction of refractive ametropia (myopia and hyperopia) and shall follow a wearing schedule as recommended by an eye care practitioner. It is also indicated to enhance color discrimination in patients with protan or duetan (red-green) color vision deficiencies. Eye Care Professional may prescribe the lens for daily wear on a frequent replacement schedule up to a year.

## Lens Insertion


1. Wash your hands thoroughly with soap and dry them with a lint free towel.
2. Place the lens on the tip of your index finger. Make sure the lens is in a bowl shape and correctly oriented.



3. Hold up your upper eyelid firmly with one finger of your other hand.



4. While holding down the lower eyelid with one finger of your applying hand, carefully place the contact lens on your cornea.



5. Slowly release upper eyelid and gently close your eyes.




6. Repeat the same procedure for the other eye.
7. When lens is fitted onto the cornea, an interval of about 20 minutes is needed for the lens to settle down.

8. In ideal condition, lens must be sitting centrally with minimal movement of 1 to 1.5mm after each blink.
9. The marking on the lens must be inserted in such a way that it is pointing at the 6 o'clock direction.

Note: Point 7 to 9 is for toric lens conditioning.

## Lens Removal

1. Wash your hands thoroughly with soap and dry them with a lint free towel.
2. Check to ensure lenses are on your cornea.
3. Hold down lower eyelid with your middle finger.
4. Use forefinger to slide lens to lower part of eye.
5. Gently pinch the bottom edge of the lens with the pads of your index finger and thumb to remove lens.



6. Repeat the same procedure for the other eye.


## Daily Cleaning and Care Regimen

1. After removing the lens from the eye, place it on the palm of your hand and apply a few drops of daily cleaner or FreshKon<sup>®</sup> Multi-Purpose Solution onto the lens.
2. Rub the lens gently with your index finger against the palm for about 20 seconds on both faces.
3. Rinse the lens thoroughly with saline or FreshKon<sup>®</sup> Multi-Purpose Solution
4. Soak the clean lens in the case with FreshKon<sup>®</sup> Multi-Purpose Solution for your next use.
5. Never rinse your lenses in water from the tap.

- Do not change the contact lens care system without consulting an Eye Care Professional.


### Schedule for Daily Wear/ Follow-up Visits

#### Schedule for Daily Wear



- Never wear your lenses beyond the period recommended by your Eye Care Professional
- Always discard lenses worn on a planned replacement schedule after the recommended wearing schedule prescribed by your Eye Care Professional. They are not indicated for overnight wear and you should not wear them while sleeping. Clinical studies have shown that the risk of serious adverse reactions is increased when these lenses are overused and worn overnight.
- This lens is meant for daily wear.





#### Follow-up visit



## Contraindications/ Warnings/ Precautions/ Adverse Reactions

1. The safety and effectiveness of contact lens depends on proper use, eye infection with possible permanent damage to vision, can result from the failure to strictly follow instructions recommended in IFU.
2. Do not use the contact lens if any of the following conditions exist: (a) Any eye disease, injury, or abnormality that affects the cornea, conjunctiva, or eyelids (b) Any active corneal infection (bacterial, fungal, protozoal, or viral) (c) If eyes becomes red or irritated.
3. Proper use and care of contact lenses and lens care products, including lens cases, are essential for the safe use of these products.
4. In the event of an adverse reaction to the wearing of the contact lens, including discomfort to the eye, red eye, excessive tearing and blurred vision and/or unusual eye secretions, the user should immediately remove the contact lens and consult an eye care professional before resuming use. It is also recommended that you see your eye care professional routinely, as directed in IFU.
5. Eye care professional should be consulted regarding the use of the contact lens in certain atmospheric or environmental conditions that can cause irritation to the eyes.
6. Avoid all harmful or irritating vapors and fumes while wearing lenses.
7. Ask your Eye Care Professional about wearing lenses during sporting activities eg. swimming.
8. Put on your lenses before putting on makeup. Do not get cosmetics, lotions, soaps, creams, deodorants, or sprays in your eyes or on your lenses.
9. Do not use if the sterile package is opened or damaged.
10. If the lens appears damaged while placing on the tip of your index finger, do not use it.
11. Do not share your contact lenses with others.
12. Always keep your lenses immersed in recommended storage solution when lenses are not being worn. Prolonged periods of drying will reduce the ability of the lens surface to return to a wetttable state.
13. Do not store the lenses or rinse the lens case with non-sterile water as this can lead to severe microbial infection. Consult Eye Care Professional for the cleaning of contact lens cases.
14. Always consult your Eye Care Professional on the suitability of lenses used and parameters.
15. The recommended storage temperature of contact lenses is below 27 °C to ensure the product shelf life is valid.

## Symbol Reference Guide on Packaging and Label

	Manufacturer	<b>CE</b> <b>1639</b>	Quality System Certification Symbol
<b>EXP</b> 	Use- by date	<b>EC</b> <b>REP</b>	Authorized representative in the European Community
<b>LOT</b>	Batch code	<b>PWR</b>	Diopter (Lens Power)
<b>STERILE</b> 	Sterilized using steam or dry heat	<b>BC</b>	Base Curve
	Consult instruction for use	<b>DIA</b>	Diameter



Printed in Singapore  
Revision date: 12/2020  
Revision number: FLEX-IFUK(3)



Obelis s.a  
Boulevard General Wahis 53  
1030 Brussels, BELGIUM  
Authorised Representative in  
the European Community



Oculus Private Limited  
10 Kaki Bukit Avenue 1  
#06-08  
Singapore 417942