

flexcon®

Type of products	Daily Wear Standard <u>Clear</u>	Daily Super Thin <u>Clear</u>	<u>Blue Tint</u> Super Thin	<u>Blue Tint UV</u> Prolong Wear	Extended Wear Standard <u>Clear</u>	Toric Blue Tint	Prosthetic Lens
Product Code	DWSD	DWST	BTST	BUPW	EWSD	TBDL	PRDW
Water content	38%	38%	38%	45%	55%	38%	38%
Material	Polymacon	Polymacon	Polymacon	Hefilcon A	Methafilcon A	Polymacon	Polymacon
Centre Thickness	0.12mm (AT -3.00D)	0.08mm (AT -3.00D)	0.08mm (AT -0.00D)	0.08mm (AT -3.00D)	0.10mm (AT -3.00D)	0.07mm (AT -0.00D)	0.12mm (AT -3.00D)
Base curve	7.80 to 9.00mm (0.10mm step)	7.80 to 9.00mm (0.10mm step)	7.80 to 9.00mm (0.10mm step)	7.80 to 9.00mm (0.10mm step)	7.80 to 9.00mm (0.10mm step)	7.80 to 9.00mm (0.10mm step)	8.40 or 8.60mm
Diameter	13.50 to 14.50mm (0.10mm step)	13.50 to 14.20mm (0.10mm step)	13.50 to 14.20mm (0.10mm step)	13.50 to 14.20mm (0.10mm step)	13.50 to 14.50mm (0.10mm step)	13.80 to 14.50mm (0.10mm step)	14.00 or 14.50mm
Power (RX)	Up to ± 20.00 (0.25D Step)	Up to ± 20.00 (0.25D Step)	Up to ± 20.00 (0.25D Step)	Up to ± 20.00 (0.25D Step)	Up to ± 20.00 (0.25D Step)	Up to ± 24.00 -1.00 to -8.00(0.25D Step) Above -8.50(0.50D Step)	Up to ± 20.00 (0.25D Step)
Axis	-	-	-	-	-	Any (In step of 5 Degree)	-